


Vascular Intervention // Coronary

Magnum

Coronary Guide Wire

Indicated for PTCA applications



- Optimal in situations with minimal distal visibility and emergency situations
- Atraumatic crossing due to olive shaped tip
- One piece core wire design provides optimal torqueability and pushability

Magnum

Coronary Guide Wire

Technical Data

Guide wire	
Diameter	0.014"
Length	175 cm
Core wire material	Stainless steel
Proximal coil	Stainless steel
Distal coil	Platinum, radiopaque
Proximal (shaft) coating	PTFE
Distal coating	Hydrophilic
Range of tip flexibilities	Between Flexible and Medium
Range of support levels	Standard
Olive tip	Diameter: 0.028" (0.7 mm) Length: 1.2 mm

Order Information

Catalogue number	Tip shape	Support level	Length
350001	Straight	Standard	175 cm
350000	J		

Magnum is part of the BIOTRONIK coronary solutions portfolio, including:

- Stents: **Orsiro, PRO-Kinetic Energy, PK Papyrus** ■ Balloons: **Pantera Lux, Pantera LE0, Pantera, AngioSculpt**
- Guide Wires: **Galeo Pro, Cruiser**

Manufacturer:
Brivant Ltd
Parkmore West Business Park
Galway
County Galway,
Ireland

Distributor:
BIOTRONIK AG
Ackerstrasse 6
8180 Bülach · Switzerland
Tel +41 (0) 44 8645111
Fax +41 (0) 44 8645005
info.vi@biotronik.com
www.biotronik.com

For ordering please contact your
local sales representative

Magnum product ordering sheet

Length 175 cm

Catalogue number	Type	Tip shape	Quantity
350001	Magnum	Straight	
350000	Magnum J	J	

Please return this sheet to:

Institute

Contact Person

Department

Address

City

Phone

Fax

E-mail

Date, Signature
