


Passeo[®]-14



Up to 3.8 x faster deflation times



Enhanced crossability



High pushability and flexibility

Passeo-14

Vascular
Intervention
Peripheral



Indicated for balloon dilatation of the stenotic portion of a lower limb artery for the purpose of improving perfusion.

Technical Data	Balloon catheter
Catheter type	OTW
Recommended guide wire	0.014"
Tip	Optimized entry profile and colored
Balloon material	SCP (Semi-Crystalline Polymer), controlled compliance (4 - 6%)
Balloon folding	3-fold
Balloon coating	Hydrophilic patchwork coating
Balloon markers	2 swaged markers (zero profile)
Sizes	ø 1.5 - 4.0 mm; L: 20 - 220 mm
Distal shaft	3.1F, hydrophilic coating, coaxial design; 150 mm length (ø 1.5/2.0 x 20 - 100 mm); 75 mm length (ø 2.0 x 140 - 220 mm and ø 2.5 - 4.0 mm)
Proximal shaft	3.9F, hydrophobic coating, coaxial design; stiffening wire
Usable length	150 cm (ø 1.5 - 4.0 mm); 120 cm (ø 1.5 - 2.0 mm); 90 cm (ø 2.5 - 4.0 mm)

Compliance Chart		Balloon diameter x length (mm)					
		ø 1.5 x 20-70	ø 2.0 x 40-220	ø 2.5 x 40-220	ø 3.0 x 40-220	ø 3.5 x 40-140	ø 4.0 x 40-140
Nominal Pressure (NP)	atm*	7	7	7	7	7	7
	ø (mm)	1.5	2.0	2.5	3.0	3.5	4.0
Rated Burst Pressure (RBP)	atm*	14	14	14	14	14	14
	ø (mm)	1.57	2.08	2.61	3.18	3.63	4.16

*1 atm = 1.013 bar

Ordering Information		Catheter Length (cm)	Balloon ø (mm)	Balloon Length (mm)						
				20	40	70	100	140	180	220
4F Antegrade approach	120	1.5	380271 ^a	380277	380283	-	-	-	-	-
	120	2.0	-	380278	380284	380290	380296	380302	380308	-
	90	2.5	-	380279	380285	380291	380297	380303	380309	-
	90	3.0	-	380280	380286	380292	380298	380304	380310	-
	90	3.5	-	380281 ^a	380287 ^a	380293 ^a	380299 ^a	-	-	-
	90	4.0	-	380282	380288	380294	380300	-	-	-
	90	4.0	-	-	-	-	-	-	-	-
4F Crossover approach	150	1.5	380313 ^a	380319	380325	-	-	-	-	-
	150	2.0	-	380320	380326	380332	380338	380344	380350	-
	150	2.5	-	380321	380327	380333	380339	380345	380351	-
	150	3.0	-	380322	380328	380334	380340	380346	380352	-
	150	3.5	-	380323 ^a	380329 ^a	380335 ^a	380341 ^a	-	-	-
	150	4.0	-	380324	380330	380336	380342	-	-	-
	150	4.0	-	-	-	-	-	-	-	-

^a8 weeks pre-order only

Passeo is a trademark or registered trademark of the BIOTRONIK Group of Companies.

BIOTRONIK AG
Ackerstrasse 6
8180 Bülach, Switzerland
Tel +41 (0) 44 8645111
Fax +41 (0) 44 8645005
info.vi@biotronik.com
www.biotronik.com

© 2019 BIOTRONIK AG – All rights reserved.
Specifications are subject to modification, revision and improvement.

BIOTRONIK
excellence for life