


Passeo[®]-18



High pushability



Low profile and wide
range of sizes



Controlled compliance

Passeo-18

Vascular
Intervention
Peripheral



Indicated to dilate stenosis in the femoral, popliteal and infrapopliteal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae.

Technical Data	Balloon catheter
Catheter type	OTW
Recommended guide wire	0.018"
Tip	Short and tapered, colored
Balloon material	SCP (Semi-Crystalline Polymer), controlled compliance (4 - 8 %)
Balloon folding	5-fold
Balloon coating	Hydrophobic patchwork coating
Balloon markers	2 swaged markers (zero profile)
Sizes	ø 2.0 - 7.0 mm; L: 20 - 200 mm
Shaft	3.8F, 3.9F (ø 6.0/7.0 mm x 170 - 200 mm); coaxial design
Usable length	90, 130 and 150 cm

Compliance Chart		Balloon diameter x length (mm)															
		ø 2.0 x 20-170	ø 2.0 x 200	ø 2.5 x 20-170	ø 2.5 x 200	ø 3.0 x 20-170	ø 3.0 x 200	ø 3.5 x 20-170	ø 3.5 x 200	ø 4.0 x 20-150	ø 4.0 x 170-200	ø 5.0 x 20-120	ø 5.0 x 150	ø 5.0 x 170-200	ø 6.0 x 20-200	ø 7.0 x 20-200	
Nominal Pressure (NP)	atm*	6	6	6	6	6	6	6	6	6	6	6	6	6	6	6	
	ø (mm)	2.0	2.0	2.5	2.5	3.0	3.0	3.5	3.5	4.0	4.0	5.0	5.0	5.0	6.0	7.0	
Rated Burst Pressure (RBP)	atm*	15	14	15	14	15	14	15	14	15	13	15	12	13	12	12	
	ø (mm)	2.1	2.1	2.6	2.6	3.2	3.2	3.7	3.7	4.3	4.2	5.3	5.2	5.2	6.2	7.3	

*1 atm = 1.013 bar

Ordering Information	Catheter Length (cm)	Balloon ø (mm)	Balloon Length (mm)								
			20	40	60	80	120	150	170	200	
			*1 atm = 1.013 bar								
Antegrade approach	4F	90	2.0	366098	366099	366100	366104	366105	366106	366114	376276
		90	2.5	357451	357458	366101	357469	357476	366107	357483	376277
		90	3.0	357452	357459	366102	357470	357477	366108	357484	376278
		90	3.5	357453	357460	366103	357471	357478	366109	357485	376279
		90	4.0	357454	357461	357465	357472	357479	366110	376272	376280
		90	5.0	357455	357462	357466	357473	357480	366111	376273	376281
		90	6.0	357456	357463	357467	357474	357481	366112	376274	376282
		90	7.0	357457	357464	357468	357475	357482	366113 ^a	376275 ^a	376283 ^a
		Retrograde approach	4F	150	2.0	366115	366118	366119	366123	366126	366129
130	2.5			357486	357491	366120	357502	357507	366130	357512	376297
130	3.0			357487	357492	366121	357503	357508	366131	357513	376298
130	3.5			357488	357493	366122	357504	357509	366132	357514	376299
130	4.0			357489	357494	357498	357505	357510	366133	376292	376300
130	5.0			357490	357495	357499	357506	357511	366134	376293	376301
130	6.0			366116	357496	357500	366124	366127	366135	376294	376302
130	7.0			366117	357497	357501	366125	366128	366136 ^a	376295 ^a	376303 ^a

^a8 weeks pre-order only

Passeo is a trademark or registered trademark of the BIOTRONIK Group of Companies.

BIOTRONIK AG
Ackerstrasse 6
8180 Bülach, Switzerland
Tel +41 (0) 44 8645111
Fax +41 (0) 44 8645005
info.vi@biotronik.com
www.biotronik.com

© 2019 BIOTRONIK AG – All rights reserved.
Specifications are subject to modification, revision and improvement.

BIOTRONIK
excellence for life